


**To obtain a  
Specimen  
Collection Kit  
contact:  
WAPHL Mailroom  
206-418-5579**

**Department of Health  
Public Health Laboratories  
1610 NE 150th Street  
Shoreline, WA 98155**



# **Chlamydia and N. gonorrhoeae Urine**

**Specimen  
Collection Kit**



**Chlamydia**

**PUBLIC HEALTH**  
ALWAYS WORKING FOR A SAFER AND  
HEALTHIER WASHINGTON

## Urine Collection Kit

## Shipping to the Public Health Laboratories

**The Gen-Probe assay** is designed to detect *C. trachomatis* and *N. gonorrhoeae* in male and female urine specimens. The PHL assay requires the use of the Gen-Probe collection kit.

Maintain specimens at 2 to 30°C and transport to the WAPHL as soon as possible..




*Transport Packages to:*

Washington Public Health  
Laboratories  
Attention: Chlamydia Lab  
1610 NE 150th Street  
Shoreline, WA 98155

### Provided in each Gen-Probe Aptima Combo 2® Collection kit:

- 50 Disposable transfer pipettes
- 50 Specimen transport tubes
- 50 Chlamydia lab forms
- Directions for use



*N. gonorrhoeae*



**Phone: 206-418-5400**

**Fax: 206-418-5445**