


**To obtain a Specimen
Collection Kit contact:
WAPHL Mailroom
206-418-5579**

**Department of Health
Public Health Laboratories
1610 NE 150th Street
Shoreline, WA 98155**



Pertussis Specimen Collection Kit



Bordetella pertussis

PUBLIC HEALTH
ALWAYS WORKING FOR A SAFER AND
HEALTHIER WASHINGTON

Pertussis Collection Kit

Instructions:

Collect the specimen with

★

Package and label according
to the Department of
Transportation regulations for
Diagnostic Specimens.

★

Place completed shipping
container into the Lab Pak
bag for USPS mailing purposes
and ship at ambient
temperature to the WAPHL
as soon as possible.

Shipping to the Public Health Laboratories



Collection Kit Contents:

- 1 Shipping container ⁶
- 1 Charcoal Media ¹
- 1 Lab Pak bag ⁷
- 2 slides ²
- 2 dacron swabs ³
- 2 nose and throat forms ⁴
- 1 pertussis immunization
history form ⁵

Transport Packages to:

Washington Public Health
Laboratories

Attention: _____
1610 NE 150th Street
Shoreline, WA 98155



Phone: 206-418-5400

Fax: 206-418-5445