


Department of Health
Public Health Laboratories
1610 NE 150th Street
Shoreline, WA 98155

**To obtain a
Specimen Collection
Kit contact:
WAPHL Mailroom
206-418-5579**



Chlamydia and N. gonorrhoeae Urethral/Cervical Specimen Collection Kit



Chlamydia


PUBLIC HEALTH
ALWAYS WORKING FOR A SAFER AND
HEALTHIER WASHINGTON

Urethral/Cervical Collection Kit

The Gen-Probe assay

is designed to detect *C. trachomatis* and *N. gonorrhoeae* in endocervical and male urethral specimens. The PHL assay requires the use of the Gen-Probe collection kit.

Maintain specimens at 2 to 30°C and transport to the WAPHL as soon as possible.



N. gonorrhoeae

Shipping to the Public Health Laboratories



Provided in each Gen-Probe Aptima Combo 2® Collection kit:

- 50 Cleaning Swabs (Female)
- 50 Unisex collection Swabs
- 50 Specimen transport tubes
- 50 Chlamydia lab forms
- Directions on Use

Transport Packages to:

Washington Public Health
Laboratories
Attention: Chlamydia Lab
1610 NE 150th Street
Shoreline, WA 98155



Phone: 206-418-5400

Fax: 206-418-5445